

MOBILE BENCHES AND CARTS



LABMODUL MOBILE BENCHES AND CARTS


 Mobility

 Storage

 Flexibility




Tray carts, depth 650 mm



Mobile benches and tray carts

In modern facilities, mobile devices can be used for safe handling. Lab module offers robust mobile tables in a wide range of applications for all purposes. Ideal for the flexible laboratory.

Mobile benches, depth 650 mm





90

labmod

Mobile benches for any purpose

The frame is supplied in powder coated steel color RAL9016 and is available with or without shelves or for under benches shelves.

The mobile benches are available in 2 standard heights or with electric height control. Height adjustment can also be supplied with battery back-up so that the benefits of the lifting function can be used anywhere without being connected to power.

The mobile table can also be delivered in stainless steel ASI316.



Mobile bench made of brushed stainless steel ASI316, with built-in electrical height control for GMO laboratories.

Low mobile benches, 2 wheels with brake (incl. Table top):

B600 x D650 x H 750 mm
B750 x D650 x H 750 mm
B900 x D650 x H 750 mm
B1200 x D650 x H 750 mm
B1500 x D650 x H 750 mm

High mobile benches, 2 wheels with brake (incl. Table top):

B600 x D650 x H 930 mm
B750 x D650 x H 930 mm
B900 x D650 x H 930 mm
B1200 x D650 x H 930 mm
B1500 x D650 x H 930 mm



Mobile bench prepared for under benches



Mobile bench prepared for under shelf



Mobile benches





Mobile benches are fitted with large, easy-running industrial wheels with high carrying capacity and stability.

Tray carts

Tray carts for 40/60 trays and baskets completes the logistics system.

Trays and baskets fit directly into cabinets and storage shelves for easy and safe handling.

The tray carts are offered in dimension: D650 x W350 or W450 x H820 mm with 4 supports

D650 x W350 or W450 x H1100 mm with 6 supports

Tower carts D650 x W450 x H1510 mm with 10 supports

Double tower carts D650 x W880 x H1510 mm with 10 supports



Mobile benches




Countertops

Material options

- Laminate
- Compact laminate
- Stainless steel
- Corian
- Glass

Laminate colors




Fronts

Materials


- Laminate
- Compact laminate
- Corian

Laminate colors




Cabinets

Low base cabinets USL



Heigh base cabinets USH




Tabletop profiles


Front edge



Postform



Compact laminate with Alu upturned edge




Trespa + laminate with Alu upturned edge



Stainless steel, with bead



Epoxy Resin with bead



Polypropylen with bead



Ceramic with bead

Rear edge

