



# labmodul

Designing safety



# Designing safety...











# HPLC CABINETS



# LABMODUL HPLC CABINETS

- ┌ Cleaner working environment
- ┌ Flexible workplace
- ┌ Integrated LED lighting
- ┌ Good insight due to Float glass
- ┌ Easy to clean
- ┌ Robust and durable





Inside the cabinet has integrated energy saving LED lighting with a pleasant light temperature.



By service and new installations, the interior cabinet goes completely out of the outer cabinet and is accessible from front, a perfect ergonomic solution.





The inner cabinet is fitted with large smooth-running industrial wheels with high stability.










## Specifications and accessories



The HPLC cabinet worktop and top shelves are provided with cable entry holes allowing an orderly and stable setup to be performed.

### Accessories

In addition to Flow Monitor, a Level Guard can be installed that continuously measures and protects against overflowing dump collectors for solvents. Level guard has built-in acoustic and visual alarm.




### Specifications

The HPLC cabinet is available in 2 sizes:  
Single unit for 1 pc. appliance  
Dimensions: WxDxH 750 x 800 x 2000 mm


Double unit for 2 pcs. appliances  
Dimensions: WxDxH 1200 x 800 x 2000 mm
















**labmodul**  
Designing safety

Version: ★ ★ ★

**LABMODUL A/S**

Valhøjs alle 190  
DK-2610 Rødovre  
[info@labmodul.dk](mailto:info@labmodul.dk)