



labmodul

Designing safety



Designing safety...









labmodul



LABMODUL SYSTEM CABINETS 60/40

- ┌ Spacious with extra storage capacity
- ┌ Prepared for trays and baskets
- ┌ High quality for daily use
- ┌ Constructed for GMO facilities
- ┌ Stylish design
- ┌ Easy to clean



System Cabinets 60/40

System cabinets are used for optimum use of space inside the storage cabinets. The cabinets can be equipped with adjustable shelves or our ISO Flex Module® baskets and trays for maximum flexibility, and efficiency.

The modular system ISO 60/40 is the result of extensive development in order to achieve maximum flexibility in systematic and safe storage that ensures order and clarity.

Universal application

System cabinets are a perfect solution for laboratories, hospitals, pharmacies, medical practices, and wherever hygienic and mobile storage solutions are needed.

A range of mobile units for the baskets and trays are available when safe transportation is needed.



A wide programme

System ISO 60/40 comes with a wide range of different material solutions for worktops and fronts and comes with different types of hinged doors or drawers, to fit exactly to the function of each cabinet.

Fixed dividers can easily be attached to the modules with a ClipOn - function and can be adjusted for every 10 and 20 mm. For a more flexible solution, the baskets and trays can be fitted with easy moveable dividers that can be placed in all positions.

Trays and baskets can be equipped with end-stops for fixing the basket or tray in the pull-out situation. When pulled out of the cabinet it can be positioned in a horizontal or sloped position. A basket or tray in a sloped position gives a good overview, even for those in the upper part of the cabinets. It is possible to add your own logo or name on trays and baskets with a separate text marker. These are offered separately on request.



Storage for DIN modules

FlexModul® stands for modular flexibility and high quality for absolutely clean storage units, for 60 cm wide cabinet modules.

Separator compartment is easily attached to the ClipOn fastening module and can be adjusted in intervals of 10 and 20 mm.

Likewise, sliding room divisions can be offered in all positions.

The trays and baskets can be equipped with end stop for fixed withdrawals. It is also possible to add its own touch on the slopes, and baskets with text or company logo.

Advantage

Separate rooms can be adapted to the content

- ◆ Integrated drawer stop
- ◆ Leakage secured
- ◆ Low weight
- ◆ ABS plastic, scratch resistant and shockproof
- ◆ Temperature resistant from -40 °C to +70 °C
- ◆ Large selection of sharing options.





Trays with over 1000 options for subdivision

Baskets and trays are available in three different materials

- 🔥 ABS plastic (gray)
- 🔥 PC plastic (transparent)
- 🔥 Wire basket (plastic coated steel)



Trays and baskets with fixed pull-out.

- 🔥 Can be placed both horizontally and sloped for easy overview of items stored.

Trays and baskets can be fitted with a 110% pull-out bearing that withstands storage up to 50 kg.

- 🔥 Ability removable rails for the absolute flexible layout of the closet



SIZES

- 600 x 400 mm
- 300 x 400 mm

HEIGHTS

- 50 mm
- 100 mm
- 200 mm

MATERIALS

PC ABS Fe

WEIGHT ALLOWANCE

50 kg

DIVISIONS

Fixed	application
Variable	application
Removable	application

Trays and baskets with pull-out can carry 20 kg to 50 kg

- 🔥 Smooth surfaces for optimum hygiene
- 🔥 Simple installation of dividers (ClipOn effect)
- 🔥 Ability labeling, example with company logo
- 🔥 Low weight
- 🔥 Stackable





HELSERUDE
SASH WINDOW

FLUJ ABRING

1/2 OF
FULL OPENING

AUTO LÄRMUNG

START/STOP


AUTO CLOSE

SYSTEM




IOI
LAB


System cabinet with ABS sides




Cabinet with solid door



Laminated door with 1 glass



Laminated door with 2 glass



Aluminum door with 2 glass





Tall storage cabinets

The individual storage system for your individual needs



Angled mount of slide possible (about 25°)



2 in 1: Straight (about 5°) and hinged (about 25°) function




Slides are easy to adjust, just push to release, max load capacity 50 kg per slide pair



Optimised for H+H FlexModules, max. storage capacity 25 H+H FlexModules à 0,186 m² = 4,65 m² per column

Angled mount of shelves with dividers and labeling possible



Alu shelves of all sizes



Trays and baskets

FlexModul® ISO formats stand for modular flexibility and high quality of absolute tidy storage in cabinet modules

Trays and baskets can, of course, be equipped with limit switches for fixed withdrawals.



See special catalog:
Mobile benches and Carts
Trays and Baskets

Materials



ABS Grey



PC transparent



Wire basket, Coated

Compartment



Fixed dividers



Variable compartments



Removable dividers
for quick cleaning

Heights



50 mm



100 mm



200 mm


Størrelser




DIN-modul
600 x 400 mm



Half DIN-modul
300 x 400 mm









labmodul
Designing safety

Version: ★ ★ ★

LABMODUL A/S

Valhøjs alle 190
DK-2610 Rødovre
info@labmodul.dk